


Discovery of diverse GPCR-targeted antibodies



The following data was presented as a poster at the US Antibody Engineering & Therapeutics Conference 2023.

INTRODUCTION

Generating diverse GPCR-targeted antibodies

There are more than 270 potential G protein-coupled receptor (GPCR) antibody drug targets spanning a wide range of indications, including oncology, autoimmunity, and metabolic disease.¹ Despite being high-value therapeutic targets, GPCRs have been largely intractable to conventional antibody discovery technologies.² Due to the deep complexity in GPCR structure and function, there are challenges at each stage of discovery.

Robust, target-specific immune responses are limited by high cross-species homology and a small number of available epitopes.² In addition, target-binding hits are extremely rare, reducing the chances of finding functional and developable leads.³

In this case study, we leveraged our discovery and development engine to find diverse antibodies against an established GPCR target, chemokine receptor 8 (CCR8). CCR8 is highly and selectively expressed on tumor-associated regulatory T cells in multiple cancer types.⁴ Antibodies against CCR8 represent a promising therapy to treat solid tumors that have progressed on standard immune checkpoint inhibitors.^{5,6,7}

To be effective, antibodies against CCR8 should be human/humanized and show favorable developability properties, highly specific binding, and potent antibody-dependent cellular cytotoxicity (ADCC) reporter activity.

Amplifying diversity, enriching for function, and optimizing developability

A combination of protein and DNA immunization campaigns was used to generate robust immune responses. We captured hundreds of diverse antibodies using deep single-cell screening. High-throughput function-based screening was used to down-select antibodies based on potency. Selected antibodies were optimized for developability using structure-guided protein engineering technologies paired with high-throughput antibody assessment.

OUTCOME

A diverse panel of anti-CCR8 antibodies optimized for developability and function

Selected molecules displayed desired properties, including:

- **Diversity:** derived from two animal species, including two humanized mouse strains
- **Potency:** strong ADCC reporter activity that is similar to a clinical benchmark
- **Developability:** favorable biophysical properties that have been de-risked by pre-emptive liability remediation

Using high-throughput immunization, protein engineering, and antibody assessment technologies, we demonstrate a strategy to identify optimized antibodies against this high-value and previously intractable class of antibody targets.

272 diverse CCR8-binding hits


To amplify antibody diversity we initiated multiple immunization campaigns in parallel, including two strains of humanized mice, two strains of rats, and diverse antigen formats (Figure 1A). We screened 2.6 million cells using a microfluidic-based single-cell screening assay and identified 633 CCR8-binding hits (Figure 1B). We progressed 272 unique sequences into high-throughput expression and characterization. This diverse set of antibodies was derived from two species and included 91 clonal families (Figure 1C).

FIGURE 1. 272 diverse CCR8-specific antibodies were identified from immunization and deep single-cell screening

(A) Two strains of humanized mice and two strains of wild-type rats were immunized using recombinant protein- or genetic-based formulations. Antibody titers were measured by flow cytometry using cells expressing human CCR8. (B) Single-cell screening assays were used to screen 2.6 million single cells and identify CCR8-specific antibodies. Representative output images of a CCR8-specific hit in a single chamber from our single-cell screen are shown. Antibodies that bind to CCR8-expressing cells were detected using fluorescence micro-scropy and images were analyzed using machine vision. (C) Antibody sequence diversity was visualized using Celium™, and 272 target-binding antibodies were progressed to high-throughput expression and characterization. Clusters represent clonal families based on V genes, J genes, and CDR3 lengths.

A Immunization strategy

Humanized mice




Rats




32
animals
screened

B Single-cell screening

Single-cell screening assay Side view



Cell screening images Top view



Legend

● Antibodies bound to CCR8-expressing cells ● Parental cells

2.6M
single
cells
screened

633
hits
recovered


C Antibody sequence diversity

Humanized mice

ATX-Gx™ 64 unique sequences | 41 clonal families
Trianni® 96 unique sequences | 18 clonal families

Wild-type rats

Sprague-Dawley 34 unique sequences | 16 clonal families
Wistar 78 unique sequences | 16 clonal families



272
unique
sequences

194 functional anti-CCR8 antibodies

We enriched the panel for antibodies with the most critical phenotype — potent ADCC reporter activity. Over 70% of expressed antibodies met the desired functional activity and we selected 24 highly potent molecules for optimization and characterization (Figure 2A).


To de-risk the panel early in development and optimize antibodies for clinical success, we identified sequence liabilities *in silico* and used this data to guide pre-emptive liability remediation (Figure 2B). The 24 selected antibodies showed a range of liabilities known to impact antibody quality, such as N-glycosylation sites, non-canonical cysteines, and deamidation motifs. In Figures 3 and 4 we describe the remediation approach of two example antibodies: (1) remediation of two sequence liabilities in a humanized mouse antibody and (2) humanization of a rat-derived antibody.

■ FIGURE 2. From 194 functional antibodies, 24 highly potent molecules were selected for high-throughput characterization and liabilities were mapped to guide pre-emptive remediation

(A) 194 of 272 expressed antibodies met the desired functional activity, measured in an ADCC reporter assay with activity quantified using luminescence. (B) Potential sequence liabilities of the 24 selected antibodies were identified *in silico* and mapped using proprietary bioinformatics pipelines. Antibody clonal diversity and liabilities were visualized using Celium™.


RLU: Relative luminescence units; FOI: Fold of induction.

(A) High-throughput screening of ADCC reporter activity



(B) *In silico* liability score overlaid with antibody sequence diversity

24 unique sequences | 17 clonal families



24 highly potent antibodies optimized for developability CONT'D


Structure-guided humanization of a highly potent rat-derived antibody

A rational approach was used to reduce predicted immunogenicity of a highly potent rat-derived antibody. We grafted complementarity-determining regions (CDRs) from the parental antibody onto a human germline framework, designed rat back-mutations in the heavy and light chains at key positions, and assessed likely impacts on structure with homology modeling (Figure 4A).

Humanized antibodies were produced and characterized in high-throughput assays to identify the optimal combination of heavy/light chain back-mutations that would reduce immunogenicity while retaining potent ADCC reporter activity. Multiple humanized variants showed reduced predicted immunogenicity, high potency, and favorable binding profiles (Figure 4B, C).

FIGURE 4. Structure-guided humanization was used to identify candidates with optimal profiles for predicted immunogenicity and function

(A) Homology modeling guided the assessment of residues targeted for humanized back-mutations on the light (L1-5) and heavy chains (H1-4). (B) The impact of heavy/light chain back-mutations on function was assessed in a rational, combinatorial approach and measured in an ADCC reporter assay as described in Figure 3B. (C) Humanization generated variants with reduced predicted immunogenicity and maintained or improved target-binding, yield, and purity compared to the parental antibody.




Anti-CCR8 antibodies with favorable developability and functional profiles that are similar to a clinical benchmark

Remediated and wild-type antibodies were expressed and characterized in high-throughput to down-select optimal antibodies. The highly diverse panel included both wild-type and remediated antibodies with favorable developability properties (Figure 5A) and potencies (Figure 5B) similar to a clinical benchmark molecule.


FIGURE 5. Strategic selection of diverse, developable antibodies with a range of potencies

(A) Selected antibodies have favorable developability profiles, including desired purity, monodispersity, relative surface hydrophobicity and heterogeneity, stability, hydrodynamics, and self-association with scores normalized to high and low controls. (B) Selected antibodies were similar to a clinical benchmark molecule in preliminary ADCC assessment. Functional activity was measured with an ADCC reporter assay as described in Figure 3B. (C) ADCC assay EC_{50} values and 95% confidence intervals are shown.

(A) Developability profiling



ADCC reporter activity



CONCLUSION

We used our antibody discovery and development engine to generate a diverse panel of anti-CCR8 antibodies that are highly potent and extensively de-risked. Our technologies and team address the unique challenges that hinder the discovery and development of antibody therapies against difficult transmembrane protein targets. Pairing a comprehensive immunization strategy with deep single B cell screening is key to identifying a large number of diverse hits against difficult GPCR targets. High-throughput developability and functional assessment, combined with structure-guided protein engineering, results in functional antibodies that are optimized for key characteristics critical to clinical success. We are leveraging these technologies against multiple GPCR and ion channel targets to unlock access to this high-value class of antibody therapeutics and treat disease areas of high unmet medical need.

REFERENCES

1. [Hutchings C, et al. \(2020\). *Antibody Therapeutics*, 3\(4\), 257-264.](#)
2. [Hutchings C, et al. \(2017\). *Nat Rev Drug Discov*, 16, 787–810.](#)
3. [Hutchings C, et al. \(2010\). *mAbs*, 2\(6\), 594-606.](#)
4. [Gallo E, et al. \(2020\). *mAbs*, 12\(1\), 1717265.](#)
5. [de Simone M, et al. \(2016\). *Immunity*, 45, 1135–1147.](#)
6. [Dépis F, et al. \(2020\). *Cancer Res*, 80 \(16 Supplement\): 4532.](#)
7. [Lu S, et al. \(2020\). *Journal for Immunotherapy of Cancer*, 8, A753–A753.](#)
8. [Lake A, et al. \(2020\). *Journal for Immunotherapy of Cancer*, 8, A769–A769.](#)