

# Breaking barriers to access intracellular targets with T-cell engagers Discovery of diverse, developable, and ultra-specific antibodies against a MAGE-A4 pMHC

#### THE OPPORTUNITY

## **Unlock intracellular tumor targets**

Immunotherapies are transforming cancer treatment by mobilizing the immune system to fight cancer. T-cell engagers — a type of immunotherapy — eliminate cancer cells by simultaneously binding tumor targets and the T-cell activating protein, CD3. Because T-cell engagers cannot access intracellular targets, their use is limited to targets expressed on the tumor cell surface.

T-cell engagers against peptides displayed on major histocompatibility complexes (pMHCs) are a promising approach for unlocking previously inaccessible, high-value, intracellular tumor targets<sup>1</sup>.

One such target is melanoma-associated antigen 4 pMHC (MAGE-A4-pMHC), a tumor-specific antigen expressed by many solid tumors, but not by most healthy tissues<sup>2</sup>.

#### THE CHALLENGE

### Finding rare antibodies against pMHCs

MAGE-A4-pMHCs are challenging immunotherapy targets because (a) proteins within the MAGE-A family are highly homologous<sup>2</sup>, (b) MAGE-A4 is often mutated in tumor cells<sup>3</sup>, and (c) there is very low expression of MAGE-A4-pMHCs on tumor cells<sup>4</sup>.

Antibodies that bind to MAGE-A4-pMHC with the specificity and affinity needed to minimize the risk of off-target binding are rare. Those that do bind and have favorable developability profiles are even more rare, making them difficult to identify using traditional approaches.

Furthermore, T-cell engager development requires tumor-binding antibodies that function as bispecifics, creating a need for diverse antibodies that can be paired with CD3-binders and tested for optimal function.

#### THE SOLUTION

## Start with diverse panels of pMHC-binding antibodies

We used our antibody discovery and development engine to discover and characterize human MAGE-A4-pMHCbinding antibodies with favorable binding and developability profiles.

Lead MAGE-A4-pMHC antibodies will be selected and paired with our panel of previously described CD3-binding antibodies using our bispecific engineering platform, OrthoMab<sup>™</sup>, to streamline the development of MAGE-A4-pMHC T-cell engagers.

#### THE OUTCOME

## We identified a panel of antibodies against MAGE-A4-pMHC with:

- diverse sequence and binding profiles
- ultra-specific, high-affinity binding
- favorable developability profiles

## 45 diverse antibodies against a MAGE-A4 peptide-MHC target

**Antibody discovery** (A) Sequence analysis



**Single-cell screening** 



Figure 1. Proprietary immunization technologies and high-throughput screening strategies were used to amplify and capture antibody diversity. (A) Humanized mice were immunized with human pMAGE-A4<sub>230-239</sub> presented on MHC-I (HLA:02\*01). (B) Multiplexed bead-based single-cell screening assays were used to screen 1.5 million single B cells and identify antibodies that bind to MAGE-A4-, but not MAGE-B2- or SARS-pMHC (an unrelated peptide control). (C) 45 MAGE-A4-pMHC-binding antibodies from 23 clonal families with sequence diversity were selected for high-throughput expression and further characterization.

# Ultra-specific, high-affinity MAGE-A4-pMHC binders

| MAGE family | Amino acid sequence | Predicted peptide affinity to MHC-I <sup>5</sup> | MHC-I relative gene expression0123 |
|-------------|---------------------|--|------------------------------------|
| A4          | GVYDG <u>REH</u> TV | SB   | -                                  |
| A10         | GLYDGMEHLI          | SB   |                                    |
| B3          | RIYDGKKHFI          | WB   | - *                                |
| B5          | QIYDGKKYYI          | WB   | *                                  |
| B6          | GIYDGILHSI          | SB   | -1                                 |
| B10         | GLYDGIEHFM          | SB   | -                                  |
| C2          | GVYAGREHFV          | WB   | *excluded due to weak              |
| sCRAP pMHC  |                     |  | peptide-WHC-I binding              |
| CELSR3      | GLSDGQWHTV          | SB   |                                    |
| LAMA1       | LLSDGKWHTV          | SB   |                                    |
| FLNC        | GLSEGHTFQV          | SB   |                                    |
| CYP4F3      | FMFEGHDTTA          | SB   | •                                  |
| VPS13B      | GLMDGSPHFL          | SB   |                                    |
| DYNC1H1     | FLSDPQVHTV          | SB   |                                    |
| FASN        | ALLDGRLQVV          | SB   |                                    |
| MMP9        | FIFEG RSYSA         | WB   |                                    |
| PIGR        |                     | SB   |                                    |
| TSPAN1      |                     | SB   | -                                  |
| TTN         | FLHDGQEYTL          | SB   |                                    |
| MALL        |                     | SB   |                                    |
| SARS        | RLQSLQTYV -         | SB   |                                    |
| non-pulsed  |                     |  | control                            |

Figure 3. MAGE-A4-pMHC-binding antibodies are highly specific. (A) Antibody binding was further assessed using a high-throughput peptide-pulsed T2 assay with MAGE-A4-pMHC and 19 other related pMHCs, including seven pMHCs from the MAGE protein family and 13 pMHCs generated using the Selective Cross-Reactive Antigen Presentation (sCRAP) algorithm<sup>6</sup>. Predicted binding affinities of peptides to MHC-I were calculated using NetMHCpan<sup>5</sup>. Relative MHC-I expression levels in T2 cells pulsed with different peptides (as a measure of positive peptide binding to MHC-I) are shown. Peptides marked with (\*) were excluded due to weak peptide-MHC-I binding.

# **Qualitative binding assessment**



Figure 2. MAGE-A4-pMHC-binding antibodies were identified. (A) Antibody binding was assessed using a multiplexed bead-based binding validation assay with MAGE-A4-pMHC and closely related pMHC antigens. Binding (> 5 FOI) was compared to a MAGE-A4-pMHC-binding benchmark antibody and detected using flow cytometry. (B) Sensorgrams generated from Surface Plasmon Resonance (SPR) kinetic analysis of a selected target-binding antibody show relatively strong binding to MAGE-A4-pMHC compared to MAGE-A8-pMHC. No binding was observed to MAGE-B2- or SARS-pMHCs.

800

#### **Cell-binding assessment**



| Pe | ptide specificity           | Cont |    |
|----|-----------------------------|------|----|
|    | pMAGE-A4 <sub>230-239</sub> | B    | В  |
|    | Non-MAGE-A4 peptides        | Μ    | Ν  |
|    |                             |      | lo |

(B) Highlighted MAGE-A4-pMHC-binding antibodies show ultra-high specificity for MAGE-A4-pMHC with low to no binding to the other pMHCs tested.



(C) Antibody binding was further assessed in a cell-based peptide-pulsed T2 assay with MAGE-A4-pMHC and closely related pMHCs. Binding (> 5 FOI) was detected using flow cytometry. Over 70% of antibodies tested bound to both MAGE-A4-pMHC and MAGE-A8-pMHC (a closely related peptide also expressed on tumor cells) in both assays. (D) Relative MHC-I gene expression levels in T2 cells pulsed with different peptides (as a measure of positive peptide binding to MHC-I) are shown.

#### AbCellera antibodies

Expressed antibodies

#### Controls

- Benchmark antibody
- Isotype control antibody

10 -

SARS

1600

Time (s

| Е-В2-рМНС | SARS-pMHC  |  |  |  |
|-----------|------------|--|--|--|
|           |            |  |  |  |
| 800 1600  | 0 800 1600 |  |  |  |

- ntrol antibodies
- Benchmark antibody MHC-binding antibody
- Isotype control antibody

# **Binding affinity** © SPR



### Antibodies

- AbCellera antibodies
- Benchmark antibody

(C) SPR kinetic analysis of antibody binding to MAGE-A4-pMHC and MAGE-A8-pMHC. No binding was observed to MAGE-B2- or SARS-pMHCs (data not shown).

# REFERENCES studies. Sci Rep-uk. 2016;6(1):25182.

- Response. Front Immunol. 2022;13:844866.
- 5. Technical University of Denmark. NetMHCpan 4.1. Pan-specific binding of peptides to MHC class I proteins of known sequence.
- https://services.healthtech.dtu.dk/services/NetMHCpan-4.1/. Accessed March 22, 2023.

#### AUTHORS Davide Tortora, Peter Bergqvist, Tim Jacobs, Grace P. Leung, Elena Vigano, Antonios Samiotakis, Harveer Dhupar, Wei Wei, Shirley Zhi, Yukiko Sato, Allison Goodman, Cindy-Lee Crichlow, Melissa Cid, Jessica Fernandes Scortecci, Ping Xiang, Ahn Lee, Kush Dalal, Vivian Li, Stephanie K. Masterman, Sherie Duncan, Aaron Yamniuk, Raffi Tonikian\*, Bryan C. Barnhart.\*presenter

DOWNLOAD POSTER 

AUTHOR AFFILIATION AbCellera, Vancouver, Canada

#1891

#### Figure 4. MAGE-A4-pMHC-binding antibodies have favorable developability profiles with high purity

measured by capillary gel electrophoresis using sodium dodecyl sulfate (CE-SDS), low aggregation measured by absolute size exclusion chromatography (aSEC), low relative surface hydrophobicity measured by analytical hydrophobic interaction chromatography (aHIC), high stability measured by nano differential scanning fluorimetry (nDSF), low polyspecificity assessed by baculovirus particle enzyme-linked immunosorbent assay (BVP-ELISA), and low self-association measured by affinity-capture self-interaction nanoparticle spectroscopy (AC-SINS) with scores normalized to high and low controls.

